Lipofectamine[®] 3000 Reagent Protocol

Protocol Outline

- A. Plate cells so they will be 70–90% confluent at the time of transfection.
- Prepare plasmid DNA-lipid complexes (recommend 2 doses of lipid).
- C. Add DNA-lipid complexes to cells.

Transfection Amounts

Component	96-well	24-well	6-well
DNA per well	100 ng	500 ng	2500 ng
P3000 [™] Reagent per well	0.2 µL	1 µL	5 µL
Lipofectamine [®] 3000 Reagent per well	0.15 and 0.3 μL	0.75 and 1.5 μL	3.75 and 7.5 μL

Transfection of siRNA

To transfect cells with siRNA, follow the protocol as described for DNA but **do not** add P3000[™] Reagent when diluting the siRNA (step 3).

Limited Product Warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale found on Life Technologies' website at www.lifetechnologies.com/termsandconditions. If you have any questions, please contact Life Technologies at www.lifetechnologies.com/support.

Important Licensing Information

These products may be covered by one or more Limited Use Label Licenses. By use of these products, you accept the terms and conditions of all applicable Limited Use Label Licenses.

Disclaimer

IFE TECHNOLOGIES CORPORATION AND/OR ITS AFFILIATE(S) DISCLAIM ALL WARRANTIES WITH RESPECT TO THIS DOCUMENT, EXPRESSED OR IMPLIED. INCLUDING BUT NOT LIMITED TO THOSE OF MERCHANTABILITY. FITNESS FOR A PARTICULAR PURPOSE. OR NON-INFRINGEMENT TO THE EXTENT ALLOWED BY LAW, IN NO EVENT SHALL LIFE TECHNOLOGIES AND/OR ITS AFFILIATE(S) BE LIABLE, WHETHER IN CONTRACT, TORT, WARRANTY, OR UNDER ANY STATUTE OR ON ANY OTHER BASIS FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT. INCLUDING BUT NOT LIMITED TO THE USE THEREOF

© 2014 Life Technologies Corporation. All rights reserved.

The trademarks mentioned herein are the property of Life Technologies Corporation and/or its affiliate(s) or their respective owners



L3000001



Publication Part No. 100022234 Publication No. MAN0009872 Rev B.0

Catalog Numbers





Size:

 $0.1 \,\mathrm{mL}$

For Research Use Only. Not for use in diagnostic procedures.

Lipofectamine® 3000 Transfection Reagent Protocol

Transfect cells according to the following table. Use the indicated volume of DNA and P3000TM Reagent with each of the two volumes of Lipofectamine® 3000 (when performing optimization). Each reaction mix volume is for one well and accounts for pipetting variations. Scale volumes proportionally for additional wells.

Timeline		eline	Steps	Procedure Details (Two Reaction Optimization)			
Day 0	4	<u>A</u>	Seed cells to be 70–90% confluent at transfection	Component	96-well	24-well	6-well
	1			Adherent cells	$1-4 \times 10^{4}$	$0.5-2 \times 10^{5}$	$0.25 - 1 \times 10^{6}$
		Diluted Lipofectamine* 3000	Dilute Lipofectamine® 3000 Reagent in Opti-MEM® Medium (2 tubes) – Mix well	Opti-MEM [®] Medium	$5 \mu L \times 2$	$25 \mu L imes 2$	$125 \ \mu L \times 2$
	2	Vortex 2-3 sec		Lipofectamine® 3000 Reagent	0.15 and 0.3 µL	0.75 and 1.5 µL	3.75 and 7.5 μL
		Diluted DNA	Prepare master mix of DNA by diluting DNA in Opti- MEM® Medium, then add P3000 [™] Reagent – Mix well	Opti-MEM [®] Medium	10 µL	50 µL	250 μL
	3			DNA (0.5–5 μg/μL)	0.2 µg	1 µg	5 µg
	J			P3000 [™] Reagent (2 µL/µg DNA)	0.4 µL	2 µL	10 µL
Day 1		Add Diluted DNA to each tube of Diluted	Diluted DNA (with P3000 [™] Reagent)	5 µL	25 μL	125 µL	
	4	4	Lipofectamine® 3000 Reagent (1:1 ratio)	Diluted Lipofectamine® 3000 Reagent	5 µL	25 µL	125 µL
	5	5	Incubate	Incubate for 5 minutes at room temperature.			
			Component (per well)	96-well	24-well	6-well	
		6	Add DNA-lipid complex to cells	DNA-lipid complex	10 µL	50 µL	250 μL
	6			DNA amount	100 ng	500 ng	2500 ng
				P3000 [™] Reagent	0.2 μL	1 µL	5 µL
				Lipofectamine® 3000 Reagent used	0.15 and 0.3 μL	0.75 and 1.5 μL	3.75 and 7.5 μL
Day 2-4	7		Visualize/analyze transfected cells	Incubate cells for 2–4 days at 37°C. Then, analyze transfected cells.			
	-				For suppor	t, visit www.lifetechn	ologies.com/support.

Invitrogen"

Lipofectamine[®] 3000—efficient, reproducible transfection for biologically-relevant cell models

Nektaria Andronikou, Xin Yu, Jia Wei, Mahalakshmi Sridharan, Rene Quintanilla, Uma Lakshmipathy, Peter Welch, Xavier de Mollerat du Jeu Life Technologies, 5791 Van Allen Way, Carlsbad, CA 92008

Introduction

With the increase of research in more biologically-relevant cell models, existing molecular and cellular techniques need to be improved. Lipofectamine® 3000, a new transfection reagent developed to improve delivery and enable use of new technologies. can be used in more relevant systems enabling faster and more reliable outcomes. Genome editing, stem cell manipulation, and immunotherapy are a few of the many rapidly growing areas that require more advanced techniques to maximize their potential applications. Lipofectamine[®] 3000 demonstrates significant improvement in a broader spectrum of cell lines when compared to current commercially available lipid-mediated transfection reagents. More importantly, Lipofectamine® 3000 has the potential to help propel many of these novel and exciting technologies forward.

Methods

Figure 1. Strategy to identify Lipofectamine[®] 3000.





Protocol

Figure 2. Enhanced and simple transfection protocol for a variety of cell lines.





Results

Figure 3. Transfection efficiency and protein expression in various cell lines.



Each reagent was used to transfect HEK 293, HeLa, LNCaP, A549, and HepG2 cells in a 96-well format, and GFP expression was analyzed 48 hours posttransfection. Lipofectamine® 3000 reagent provided higher GFP transfection efficiency than Lipofectamine® 2000 and FuGENE® HD reagents for all five cell lines.

Figure 4. Delivery of siRNA for gene knockdown.



Lipofectamine[®] 3000 is a versatile reagent that can also be used to deliver siRNA using the same transfection protocol. Simply substitute siRNA for DNA. For this experiment, knockdown of endogenous luciferase was achieved in three engineered luciferase cell lines using Lipofectamine[®] 3000. Lipofectamine® 2000, and Lipofectamine® RNAiMax. Reagents were complexed with Silencer® Select siRNA targeting luciferase at the specified siRNA dosages.

Figure 5. Transfection in H9 embryonic stem cells.



Transfection performed in H9 human embryonic stem cells with Lipofectamine® 2000 and Lipofectamine® 3000 in a 96-well format. GFP expression analysis was performed 24 hours posttransfection.

Table 1. Culture and transfection conditions Figure 6. Reprogramming results

Delivery	Cell	Media change
method	density	posttransfection
Neon [®]	1.0 x 10 ⁴	Day 1–Day 14: use
Transfection	per 6-well	N2B27 media with
System	plate	100 ng/mL FGF
Lipofectamine® 3000	3.0 x 104 per 6-well plate	Day 15–Day 16: use StemPro® SFM Day 17: stain with alkaline phosphatase

Cell culture and transfection conditions. Transfection was performed in BJ fibroblasts using the Neon® Transfection System at the recommended conditions and Lipofectamine[®] 3000 at 3.6 µL per well. Epi5™ Episomal iPSC Reprogramming Vectors were used (Cat. No. A14703). Media changes were performed daily according to the following protocol: Generation of human induced pluripotent stem cells (hiPSCs) from fibroblasts using episomal vectors (see Table 1).

Figure 7. Genomic modification using GeneArt® CRISPR Nuclease Vector Kits



The all-in-one GeneArt® CRISPR vector system contains a Cas9 nuclease expression cassette and a guide RNA cloning cassette that was used to target the AAVS1 safe harbor locus; a downstream orange fluorescent protein (OFP) gene helps determine delivery efficiency and can also be used for enrichment. Lipofectamine[®] 2000 and Lipofectamine[®] 3000 were used to transfect U2OS and HepG2 cells in a 12-well format. Efficiency and OFP expression were analyzed 72 hours posttransfection and (A) U2OS and (B) HepG2 cells showed 4-fold and 80-fold improvement with Lipofectamine® 3000, respectively.

Α



Cleavage efficiency determined with the GeneArt® Genomic Cleavage Detection Kit. Lipofectamine® 2000 and Lipofectamine® 3000 were used to deliver GeneArt® Precision TALs and CRISPR nucleases targeting the AAVS1 safe harbor locus in U2OS and HepG2 cell lines in a 12-well format. Cell lysates were collected and processed to determine cleavage efficiency. Increased TALEN- and CRISPR-mediated cleavage were observed in both cell lines transfected with Lipofectamine® 3000. (A) U2OS cells transfected with Lipofectamine® 3000 showed 1.5-fold improved TALEN cleavage efficiency and slightly improved CRISPR cleavage. (B) HepG2 cells had 3-fold and 8-fold improved efficiency for TALEN- and CRISPR-mediated cleavage, respectively.



Lipofectamine[®] 3000 reagent for iPSC reprogramming.



Results obtained, via brightfield microscopy, for (A) Lipofectamine® 3000 and (B) the Neon® Transfection System indicate that reprogramming was successful in generating iPSC colonies. A terminal stain was performed with red alkaline phosphatase and colony counts are indicated in (C) and (D). Detailed culture protocols can be found at lifetechnoloiges.com



Figure 8. Genomic cleavage detection of the AAVS1 safe harbor locus.

Cell type	Lipofectamine® 3000 reagent transfection efficiency	Fold protein expression improvement, Lipofectamine® 3000 vs. 2000 reagent	Cell type	Lipofectamine® 3000 reagent transfection efficiency	Fold protein expression improvement, Lipofectamine® 3000 vs. 2000 reagent
3Т3		4	L6		8
4T1		2	L929		2
A431		2	LNCaP		6
A549		3	MCF 10A		5
ACHN		2	MCF7		2
bEnd.3		9	MDA-MB-231		3
BJ		3	MDA-MB-435		1
BT-549		4	MDA-MB-468		9
C2C12		3	MDCK		1
C6		5	Neuro-2a		1
Caco-2		2	NCI-H23		2
Caki-1		4	NCI-H460		17
СНО-К1		1	P19		1
CHO-S		1	PANC-1		3
COLO 205		4	PC12		2
COS-7		4	RAW264.7		4
DU 145		2	RBL-2H3		2
H460		3	RD		4
H9c2		3	Saos-2		4
HCC1937		5	SH-SY5Y		1
HCT116		1	SK-BR-3		4
HEK 293		2	SK-MEL-28		2
HeLa		3	SK-N-SH		6
Hep-3B		2	SK-OV-3		3
Нера 1-6		1	SW480		2
HepG2		9	SW620		5
Hs 578T		3	T98G		4
cHT29		1	U20S		3
Huh-7		4	U937		2
Jurkat		1	Vero		1
K-562		1			
Transfection efficiency (%): <30	% 30-50% 51-7	9% >80%			

